







Tube #	Salt	Tube #	Polymer	Tube #	pH \diamond
A1	1. 0.2 M Sodium fluoride	1.	20% w/v Polyethylene glycol 3,350	1.	7.3
A2	2. 0.2 M Potassium fluoride	2.	20% w/v Polyethylene glycol 3,350	2.	7.3
A3	3. 0.2 M Ammonium fluoride	3.	20% w/v Polyethylene glycol 3,350	3.	6.2
A4	4. 0.2 M Lithium chloride	4.	20% w/v Polyethylene glycol 3,350	4.	6.8
A5	5. 0.2 M Magnesium chloride hexahydrate	5.	20% w/v Polyethylene glycol 3,350	5.	5.9
A6	6. 0.2 M Sodium chloride	6.	20% w/v Polyethylene glycol 3,350	6.	6.9
A7	7. 0.2 M Calcium chloride dihydrate	7.	20% w/v Polyethylene glycol 3,350	7.	5.1
A8	8. 0.2 M Potassium chloride	8.	20% w/v Polyethylene glycol 3,350	8.	7.0
A9	9. 0.2 M Ammonium chloride	9.	20% w/v Polyethylene glycol 3,350	9.	6.3
A10	10. 0.2 M Sodium iodide	10.	20% w/v Polyethylene glycol 3,350	10.	7.0
A11	11. 0.2 M Potassium iodide	11.	20% w/v Polyethylene glycol 3,350	11.	7.0
A12	12. 0.2 M Ammonium iodide	12.	20% w/v Polyethylene glycol 3,350	12.	6.2
B1	13. 0.2 M Sodium thiocyanate	13.	20% w/v Polyethylene glycol 3,350	13.	6.9
B2	14. 0.2 M Potassium thiocyanate	14.	20% w/v Polyethylene glycol 3,350	14.	7.0
B3	15. 0.2 M Lithium nitrate	15.	20% w/v Polyethylene glycol 3,350	15.	7.1
B4	16. 0.2 M Magnesium nitrate hexahydrate	16.	20% w/v Polyethylene glycol 3,350	16.	5.9
B5	17. 0.2 M Sodium nitrate	17.	20% w/v Polyethylene glycol 3,350	17.	6.8
B6	18. 0.2 M Potassium nitrate	18.	20% w/v Polyethylene glycol 3,350	18.	6.8
B7	19. 0.2 M Ammonium nitrate	19.	20% w/v Polyethylene glycol 3,350	19.	6.2
B8	20. 0.2 M Magnesium formate dihydrate	20.	20% w/v Polyethylene glycol 3,350	20.	7.0
B9	21. 0.2 M Sodium formate	21.	20% w/v Polyethylene glycol 3,350	21.	7.2
B10	22. 0.2 M Potassium formate	22.	20% w/v Polyethylene glycol 3,350	22.	7.3
B11	23. 0.2 M Ammonium formate	23.	20% w/v Polyethylene glycol 3,350	23.	6.6
B12	24. 0.2 M Lithium acetate dihydrate	24.	20% w/v Polyethylene glycol 3,350	24.	7.9
C1	25. 0.2 M Magnesium acetate tetrahydrate	25.	20% w/v Polyethylene glycol 3,350	25.	7.9
C2	26. 0.2 M Zinc acetate dihydrate	26.	20% w/v Polyethylene glycol 3,350	26.	6.4
C3	27. 0.2 M Sodium acetate trihydrate	27.	20% w/v Polyethylene glycol 3,350	27.	8.0
C4	28. 0.2 M Calcium acetate hydrate	28.	20% w/v Polyethylene glycol 3,350	28.	7.5
C5	29. 0.2 M Potassium acetate	29.	20% w/v Polyethylene glycol 3,350	29.	8.1
C6	30. 0.2 M Ammonium acetate	30.	20% w/v Polyethylene glycol 3,350	30.	7.1
C7	31. 0.2 M Lithium sulfate monohydrate	31.	20% w/v Polyethylene glycol 3,350	31.	6.0
C8	32. 0.2 M Magnesium sulfate heptahydrate	32.	20% w/v Polyethylene glycol 3,350	32.	6.0
C9	33. 0.2 M Sodium sulfate decahydrate	33.	20% w/v Polyethylene glycol 3,350	33.	6.7
C10	34. 0.2 M Potassium sulfate	34.	20% w/v Polyethylene glycol 3,350	34.	6.8
C11	35. 0.2 M Ammonium sulfate	35.	20% w/v Polyethylene glycol 3,350	35.	6.0
C12	36. 0.2 M Sodium tartrate dibasic dihydrate	36.	20% w/v Polyethylene glycol 3,350	36.	7.3
D1	37. 0.2 M Potassium sodium tartrate tetrahydrate	37.	20% w/v Polyethylene glycol 3,350	37.	7.4
D2	38. 0.2 M Ammonium tartrate dibasic	38.	20% w/v Polyethylene glycol 3,350	38.	6.6
D3	39. 0.2 M Sodium phosphate monobasic monohydrate	39.	20% w/v Polyethylene glycol 3,350	39.	4.7
D4	40. 0.2 M Sodium phosphate dibasic dihydrate	40.	20% w/v Polyethylene glycol 3,350	40.	9.1
D5	41. 0.2 M Potassium phosphate monobasic	41.	20% w/v Polyethylene glycol 3,350	41.	4.8
D6	42. 0.2 M Potassium phosphate dibasic	42.	20% w/v Polyethylene glycol 3,350	42.	9.2
D7	43. 0.2 M Ammonium phosphate monobasic	43.	20% w/v Polyethylene glycol 3,350	43.	4.6
D8	44. 0.2 M Ammonium phosphate dibasic	44.	20% w/v Polyethylene glycol 3,350	44.	8.0
D9	45. 0.2 M Lithium citrate tribasic tetrahydrate	45.	20% w/v Polyethylene glycol 3,350	45.	8.4
D10	46. 0.2 M Sodium citrate tribasic dihydrate	46.	20% w/v Polyethylene glycol 3,350	46.	8.3
D11	47. 0.2 M Potassium citrate tribasic monohydrate	47.	20% w/v Polyethylene glycol 3,350	47.	8.3
D12	48. 0.2 M Ammonium citrate dibasic	48.	20% w/v Polyethylene glycol 3,350	48.	5.1

\diamond Measured pH at 25 °C



PEG/Ion Screen contains forty-eight unique reagents. To determine the formulation of each reagent, simply read across the page.

HAMPTON
RESEARCH

Solutions for Crystal Growth

© 1991-2010 Hampton Research Corp. all rights reserved
Printed in the United States of America. This guide or
parts thereof may not be reproduced in any form without
the written permission of the publishers.

Tube #	Salt	Tube #	Buffer ◇	Tube #	Polymer
E1	1. 0.1 M Sodium malonate pH 4.0	1.	None	1.	12% w/v Polyethylene glycol 3,350
E2	2. 0.2 M Sodium malonate pH 4.0	2.	None	2.	20% w/v Polyethylene glycol 3,350
E3	3. 0.1 M Sodium malonate pH 5.0	3.	None	3.	12% w/v Polyethylene glycol 3,350
E4	4. 0.2 M Sodium malonate pH 5.0	4.	None	4.	20% w/v Polyethylene glycol 3,350
E5	5. 0.1 M Sodium malonate pH 6.0	5.	None	5.	12% w/v Polyethylene glycol 3,350
E6	6. 0.2 M Sodium malonate pH 6.0	6.	None	6.	20% w/v Polyethylene glycol 3,350
E7	7. 0.1 M Sodium malonate pH 7.0	7.	None	7.	12% w/v Polyethylene glycol 3,350
E8	8. 0.2 M Sodium malonate pH 7.0	8.	None	8.	20% w/v Polyethylene glycol 3,350
E9	9. 4% v/v Tacsimate pH 4.0	9.	None	9.	12% w/v Polyethylene glycol 3,350
E10	10. 8% v/v Tacsimate pH 4.0	10.	None	10.	20% w/v Polyethylene glycol 3,350
E11	11. 4% v/v Tacsimate pH 5.0	11.	None	11.	12% w/v Polyethylene glycol 3,350
E12	12. 8% v/v Tacsimate pH 5.0	12.	None	12.	20% w/v Polyethylene glycol 3,350
F1	13. 4% v/v Tacsimate pH 6.0	13.	None	13.	12% w/v Polyethylene glycol 3,350
F2	14. 8% v/v Tacsimate pH 6.0	14.	None	14.	20% w/v Polyethylene glycol 3,350
F3	15. 4% v/v Tacsimate pH 7.0	15.	None	15.	12% w/v Polyethylene glycol 3,350
F4	16. 8% v/v Tacsimate pH 7.0	16.	None	16.	20% w/v Polyethylene glycol 3,350
F5	17. 4% v/v Tacsimate pH 8.0	17.	None	17.	12% w/v Polyethylene glycol 3,350
F6	18. 8% v/v Tacsimate pH 8.0	18.	None	18.	20% w/v Polyethylene glycol 3,350
F7	19. 0.1 M Succinic acid pH 7.0	19.	None	19.	12% w/v Polyethylene glycol 3,350
F8	20. 0.2 M Succinic acid pH 7.0	20.	None	20.	20% w/v Polyethylene glycol 3,350
F9	21. 0.1 M Ammonium citrate tribasic pH 7.0	21.	None	21.	12% w/v Polyethylene glycol 3,350
F10	22. 0.2 M Ammonium citrate tribasic pH 7.0	22.	None	22.	20% w/v Polyethylene glycol 3,350
F11	23. 0.1 M DL-Malic acid pH 7.0	23.	None	23.	12% w/v Polyethylene glycol 3,350
F12	24. 0.2 M DL-Malic acid pH 7.0	24.	None	24.	20% w/v Polyethylene glycol 3,350
G1	25. 0.1 M Sodium acetate trihydrate pH 7.0	25.	None	25.	12% w/v Polyethylene glycol 3,350
G2	26. 0.2 M Sodium acetate trihydrate pH 7.0	26.	None	26.	20% w/v Polyethylene glycol 3,350
G3	27. 0.1 M Sodium formate pH 7.0	27.	None	27.	12% w/v Polyethylene glycol 3,350
G4	28. 0.2 M Sodium formate pH 7.0	28.	None	28.	20% w/v Polyethylene glycol 3,350
G5	29. 0.1 M Ammonium tartrate dibasic pH 7.0	29.	None	29.	12% w/v Polyethylene glycol 3,350
G6	30. 0.2 M Ammonium tartrate dibasic pH 7.0	30.	None	30.	20% w/v Polyethylene glycol 3,350
G7	31. 2% v/v Tacsimate pH 4.0	31.	0.1 M Sodium acetate trihydrate pH 4.6	31.	16% w/v Polyethylene glycol 3,350
G8	32. 2% v/v Tacsimate pH 5.0	32.	0.1 M Sodium citrate tribasic dihydrate pH 5.6	32.	16% w/v Polyethylene glycol 3,350
G9	33. 2% v/v Tacsimate pH 6.0	33.	0.1 M BIS-TRIS pH 6.5	33.	20% w/v Polyethylene glycol 3,350
G10	34. 2% v/v Tacsimate pH 7.0	34.	0.1 M HEPES pH 7.5	34.	20% w/v Polyethylene glycol 3,350
G11	35. 2% v/v Tacsimate pH 8.0	35.	0.1 M Tris pH 8.5	35.	16% w/v Polyethylene glycol 3,350
G12	36. None	36.	0.07 M Citric acid, 0.03 M BIS-TRIS propane / pH 3.4	36.	16% w/v Polyethylene glycol 3,350
H1	37. None	37.	0.06 M Citric acid, 0.04 M BIS-TRIS propane / pH 4.1	37.	16% w/v Polyethylene glycol 3,350
H2	38. None	38.	0.05 M Citric acid, 0.05 M BIS-TRIS propane / pH 5.0	38.	16% w/v Polyethylene glycol 3,350
H3	39. None	39.	0.04 M Citric acid, 0.06 M BIS-TRIS propane / pH 6.4	39.	20% w/v Polyethylene glycol 3,350
H4	40. None	40.	0.03 M Citric acid, 0.07 M BIS-TRIS propane / pH 7.6	40.	20% w/v Polyethylene glycol 3,350
H5	41. None	41.	0.02 M Citric acid, 0.08 M BIS-TRIS propane / pH 8.8	41.	16% w/v Polyethylene glycol 3,350
H6	42. 0.02 M Calcium chloride dihydrate, 0.02 M Cadmium chloride hydrate, 0.02 M Cobalt(II) chloride hexahydrate	42.	None	42.	20% w/v Polyethylene glycol 3,350
H7	43. 0.01 M Magnesium chloride hexahydrate 0.005 M Nickel(II) chloride hexahydrate	43.	0.1 M HEPES sodium pH 7.0	43.	15% w/v Polyethylene glycol 3,350
H8	44. 0.02 M Zinc chloride	44.	None	44.	20% w/v Polyethylene glycol 3,350
H9	45. 0.15 M Cesium chloride	45.	None	45.	15% w/v Polyethylene glycol 3,350
H10	46. 0.2 M Sodium bromide	46.	None	46.	20% w/v Polyethylene glycol 3,350
H11	47. 1% w/v Tryptone	47.	0.05 M HEPES sodium pH 7.0	47.	12% w/v Polyethylene glycol 3,350
H12	48. 1% w/v Tryptone	48.	0.05 M HEPES sodium pH 7.0	48.	20% w/v Polyethylene glycol 3,350

◇ Buffer pH is that of a 1.0 M stock prior to dilution
with other reagent components: pH with HCl or NaOH.

PEG/Ion 2 Screen contains forty-eight unique reagents. To determine the formulation of each reagent, simply read across the page.

HAMPTON
RESEARCH

Solutions for Crystal Growth